

MICROCOPY RESOLUTION TEST CHART NATIONAL BUREAU OF STANDARDS-1963-A





ADE00050

NRL Memorandum Report 5352

Software for a Computer Controlled Laser Doppler Velocimeter

C. R. KAPLAN, S. R. LUSTIG AND F. W. WILLIAMS

Combustion and Fuels Branch
Chemistry Division

June 29, 1984







NAVAL RESEARCH LABORATORY Washington, D.C.

Approved for public release; distribution unlimited.

84 07 27 081




REPORT DOCUMENTATION PAGE							
14 REPORT SECURITY CLASSIFICATION UNCLASSIFIED		TE RESTRICTIVE MARKINGS					
2& SECURITY CLASSIFICATION AUTHORITY		3. DISTRIBUTION/AVAILABILITY OF REPORT					
26 DECLASSIFICATION/DOWNGRADING SCHEDULE		Approved for public release; distribution unlimited.					
4. PERFORMING ORGANIZATION REPORT NUMBER(S)		S. MONITORING ORGANIZATION REPORT NUMBERIS)					
NRL Memorandum Report 5352							
64 NAME OF PERFORMING ORGANIZATION	Bb. OFFICE SYMBOL 11f applicable:	7a. NAME OF MONITORING ORGANIZATION					
Naval Research Laboratory	Code 6180						
6c ADDRESS (City State and ATP Code)		7b. ADDRESS (City, State and ZIP Code)					
Washington, DC 20375							
8. NAME OF FUNDING/SPONSORING ORGANIZATION	86 OFFICE SYMBOL (If applicable)	9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER					
Naval Material Command							
Sc. ADDRESS (City, State and ZIP Code)		10. SOURCE OF FU	NDING NOS.				
Washington, DC 20360		PROGRAM ELEMENT NO.	PROJECT	TASK NO.	WORK UNIT NO.		
11 TITLE Include Security Classification: Softw	ana fan	62543N	SF43-400-		DN120-014		
Computer Controlled Laser Doppler			391				
12. PERSONAL AUTHOR(S)							
Kaplan, C.R., Lustig, S.R., and Willi		14. DATE OF REPO	BT /Ve Me Dev	15 PAGE CO	NINT.		
	5/82 to 8/82	June 29, 19	-	32	,0		
16. SUPPLEMENTARY NOTATION							
17 COSATI CODES	18. SUBJECT TERMS (C	ontinue on reverse if n	ecessary and identi	fy by block numbers			
FIELD GROUP SUB. GR.	Laser doppler ve			Combustion			
	Turbulence						
Software, consisting of a main FORTRAN program and an assembly language driver, has been developed for the computer control of a laser doppler velocimeter (LDV). This is part of a continuing effort in the use of the LDV to measure turbulence in combustion processes. A brief discussion of the LDV hardware is included. An illustrative example follows to demonstrate software capabilities.							
20 DISTRIBUTION/AVAILABILITY OF ABSTRACT		21 ABSTRACT SECURITY CLASSIFICATION UNCLASSIFIED					
UNCLASSIFIED UNLIMITED 🖾 SAME AS RPT 🗀 DTIC USERS 🗆		226 TELEPHONE N		22c OFFICE SYME	101		
228 NAME OF RESPONSIBLE INDIVIDUAL		Anciede Area Co	ode:				
C. R. Kaplan		(202) 767-2	2476	Code 618	0		




CONTENTS

I.	INTRODUCTION	. 1
II.	EQUIPMENT	. 2
	A. Laser B. Counter Processor C. Buffer Interface	. 2
ш.	SOFTWARE	. 5
	A. Driver B. Program LDV	
IV.	EXAMPLE RUN	. 8
v.	CONCLUSIONS	. 9
VI.	REFERENCES	. 9
APP	PENDIX A — SOFTWARE DRIVER	. 11
APP	PENDIX B — MAIN FORTRAN PROGRAM	. 14
APP	PENDIX C — EXAMPLE: SIX cm METHANOL PAN FIRE	. 25







SOFTWARE FOR A COMPUTER CONTROLLED LASER DOPPLER VELOCIMETER

I. INTRODUCTION

The hazard associated with fire in confined spaces is critical to the Navy due to the vulnerability of a ship to the disastrous effects of an uncontrollable fire. An understanding of the physical mixing processes in diffusion flames, which are characteristic of most unwanted fires, is necessary. By knowing the mixing patterns, a better understanding of fire growth in confined areas can be developed. This ultimately will lead to ways of controlling materials, geometry, design and structure to fire harden our ships more. The need for better mixing information has prompted research efforts toward the development of a technique to measure turbulence in combustion processes using laser doppler velocimetry (LDV).

The LDV is a state-of-the-art diagnostic tool to measure the velocity of particles moving in a flow field. The ability to characterize highly turbulent flow fluctuations at near zero mean velocity and to make measurements without obstructing the flow makes turbulence measurement via LDV preferable to that of other developed methods, such as the hot wire anemometer. The inability of the hot wire anemometer to measure turbulence in smoke-laden environments makes the LDV approach particularly attractive (1,2).

The principle behind laser doppler velocimetry lies in the phenomenon that particles suspended in a flow field will scatter laser light that is directed upon them. The frequency of the scattered light depends on the particle velocity and the angle through which the light has been scattered. The frequency of the scattered light is different from the frequency of the incident laser light. This frequency difference is the doppler frequency, which can be detected by electronic processing.

This is a continuing research effort, the result of which, to date, has been the development of computer software to control the LDV system. An example run follows, for illustrative purposes in demonstration of the software capabilities.

Manuscript approved March 22, 1984.

II. EQUIPMENT

The LDV System (Disa Electronics) is comprised of a two-color laser and associated optics, two counter processors, electronic frequency shifters, and a buffer interface. Figure 1 illustrates the overall schematic of the LDV system. The pieces of the equipment describing the LDV System are described below.


II. A. LASER

The blue/green colored beam from a four watt, argon ion laser is passed through a series of color-neutral and color-separating beam splitters to create three beams, blue/green, green and blue, as shown in Figure 2. The three beams are focused by the transmitting optics to an intersection point, approximately 16 feet from the laser where the fluid whose flow patterns are being measured is located. The purpose of choosing such an unusually large focal length (most LDV systems employ a three to four foot focal length) is to have the capability, when necessary, to make measurements during large scale fire tests where it is desirable to keep the hardware at a sufficiently large distance from the fire source.


The intersection of the green and blue/green beams results in the formation of interference fringes. A particle from the flow field which passes through the interference fringes scatters light, the intensity of which rises and falls as the particle moves from bright fringe to dark fringe. scattered light is collected on a photomultiplier tube in the backscatter mode and converted to a frequency-modulated current burst. Further processing occurs in the counter processors (as discussed below) to give the vertical component of velocity of the particle in the flow field. The intersection of the blue and blue/green beams similarly creates interference fringes from which a particle passing through scatters This scattered light is collected on a second photomultiplier tube, converted to a doppler signal, and further processed in the second counter processor to give the horizontal velocity component of the same particle in the flow field.

II. B. COUNTER PROCESSOR

The frequency-modulated current burst from each of the photomultipliers is amplified and band-pass filtered in separate counter processors. Each burst represents the passage of a single particle and contains information on the number of fringes passed by the particle. The frequency of the doppler burst signal is detected by internal logic circuitry. The distance between the fringes is known from the geometry of the optical beams and the wavelength of the laser



Schematic of LDV System Figure 1.



separate blue, green and combination blue/green laser Optical arrangement in backscatter mode resulting in beams, which intersect to form interference fringes. Particles in flow-field scatter light, the intensity of which rises and falls as the particle moves from bright fringe to dark fringe. Figure 2.

line. Instantaneous particle velocity may then be calculated by multiplying the fringe spacing by the measured doppler signal frequency.

The counter processors have the capability to send four types of data to the buffer interface. Doppler frequency data and sample interval time (time interval between each successive validated signal) data are contained on the pins of one of the rear panel multi-connectors of the counter processor. The other gives information on the number of fringes and burst time, i.e., time between each doppler burst.

II. C. BUFFER INTERFACE

The Buffer Interface acts as a data and control interface in the transfer of digital data between the counter processors and the Hewlett-Packard-1000 mini-computer. It functions as a buffer memory for the randomly arriving data from the counter processors before it is transferred to the HP-1000. It also receives the control words from the computer, as written in the software, so that the LDV will operate with the appropriate data and control channels opened.

Synchronization within a selected time interval between the data from the two counter processors is achieved with the coincidence filter board inside the Buffer Interface. The selection of input and output channels and the setting of coincidence time intervals are all software-controlled.

The Buffer Interface has been designed to interface with the DEC PDP-11 mini-computer family. An additional circuit board was added to the Buffer Interface which converts the control and data format to that of the HP Microcircuit Board 12566B for Hewlett Packard mini-computer in the 2100 series.

Data are transferred from the additional interface board, with a 16 bit, parallel full duplex format, to the 12566 Board of the HP 21MX series mini-computer. The 12566 Board has 16 bit input and output registers. A system of electronic "handshakes" controls the signal and data lines, as described in the software section of this report.

III. SOFTWARE

The transfer of information to and from the LDV is software controlled. Control of the LDV by the HP-1000 computer is achieved by a software driver. The driver is called by the main FORTRAN program, program LDV, both of which were written at the Naval Research Laboratory (NRL) and are discussed in the following sections.

III. A. DRIVER

A relocatable HP-assembly language driver was written to enable the computer to communicate with the LDV Buffer Interface. The driver, as shown in Appendix A, is FORTRAN callable, and contains all of the necessary "set control," "set flag," etc. instructions to control and read data from the LDV.

A non-interrupt data transfer method was chosen due to the higher degree of ease with which this type of driver could be written compared to that of an interrupt-type driver. non-interrupt method involves a "wait-for-flag" routine when inputting LDV data into the data-register. The computer signals the interface card to start the LDV device. driver loops until the LDV is ready to input data into the data-register. Readiness to input data is expressed by setting a flag flip-flop. After the flag is set, the program continues. At times, the doppler signal is of insufficient strength to cause the flag to set. To prevent the driver from looping indefinitely, a routine to wait 1000 times for the flag was incorporated. If the flag does not set after 1000 attempts, the program is terminated. The wait-for-flag routine is not required when outputting control words from the computer to the LDV.

The memory protect logic feature, which is incorporated in the HP-1000 Real Time Executive operating system, limits control of all Input/Output operations to interrupt control only. Thus, the memory protect fence is disabled by assembly language software, so that the non-interrupt driver has control of the LDV.

The driver consists of two main segments, the initiation section and the continuation section. The initiation section performs two basic functions: (a) to clear the buffer memory before a new record of data is to be entered, and (b) to send the control words to the LDV.

The collection of data is accomplished in the continuation section. Counting variables are incremented each time a data word is collected so that data rate can be acknowledged. The memory protect fence is re-enabled after a buffer of data is collected. This gives the interrupt system control of the input/output peripheral devices before returning to main FORTRAN program. Three control words are sent to the LDV which open up appropriate channels for data collection. The control words, which are collected in the initiation section of the driver and are sent out in the main program, perform three functions.

The first word turns on the coincidence filter board and opens up channels 1 and 2. The second control word sets the time interval between successive samples to 125 microseconds. This will ensure that the data from the two counter processors will be synchronized within this time interval. The third control word is the static control word and must remain set during operation. Channels 1 and 2 (connected with the counter processors) and channel 4 (connected with the Coincidence Filter Board) are opened for input. Doppler frequency data and sample interval data may be simultaneously collected from one of the multi-connectors on the counter processor; fringe number and burst time data are simultaneously collected from the other.

Because the 12566 Board does not give a data transmitted signal after a data word has been read into the computer, this signal is created artificially by transmitting a device command signal with the address of octal 60000. In both the initiation and continuation sections, the number 60000 (octal) is loaded into the B register. The instructions OTB LDV (output the contents of the B register to the LDV) and STC LDV (set the control flip flop on the LDV interface card in the computer I/O cage) result in the artificial creation of a data transmitted signal.

III. B. PROGRAM LDV

The main FORTRAN program, as shown in Appendix B, is comprised of three main sections which contain calls to the assembly language driver and calls to several subroutines. The first segment is used to test the LDV by taking samples and writing this out to any terminal. The user may specify any control words that are desired. This is a good way to test to see if the different channels which are opened are actually sending data.

The second segment is used to collect data for a userspecified amount of time and writes the raw data to a created
disc file. The maximum number of 512 data transfer words can
be collected in a record, as this is the size of buffer memory
in the Buffer Interface. The number of data words collected,
represented by the variable ICNT, is one of the variables
passed to the driver from the main program, which enables the
calculation of data rate.

Segment number three is used to sort and assemble the raw data collected on the disc file created by segment number two. A second data file is created which contains the sorted data. This segment reads in a record of data of up to 512 words from the raw data file. Subroutine ASMBL, which is then called, evaluates each data word based on the type of information that corresponds to each of the 16 bits comprising the data words.

This information is passed to subroutine SORT which identifies the data as either doppler frequency data, sample interval data, number of fringes data, or burst time data. Doppler frequency data is converted to velocity, by the following relationship:

Velocity = Doppler Frequency x Fringe Spacing

where

Fringe Spacing =
$$\frac{\lambda}{2 \sin \theta/2}$$

 λ = laser wavelength

 θ = angle of intersection of laser beams

Summations of the velocity data are performed in Subroutine STAT for the purpose of calculating average velocity at a subsequent point in the program. Subroutine GATHR groups the velocity and sample interval time data by burst and writes the assembled data to the output file. At this point, program statement number 310, the program loops back up to program statement number 309 to read the next record of data of up to 512 words from the raw data file.

After the final record of data has been written to the output file, Subroutine STAT calculates average velocity for each of the counters, i.e., horizontal and vertical component, and the magnitude and direction of the resultant velocity vector.

IV. EXAMPLE RUN

The following example is intended as an illustration of the utility of the software in transferring information to and from the LDV.

The vertical and horizontal components of velocity were measured in a test run for one point in a 6 cm methanol pan fire. The laser beam intersected at a vertical position 1 cm above the lip of the pan and in the center of the pan in the horizontal direction. The liquid level of methanol was maintained constant (3). Data were collected for approximately 30 minutes.

The lack of formation of soot particles in the methanol flame, due to methanol's clean burning characteristics, creates some difficulty with LDV measurement. The principle behind the LDV technique rests in the fact that particles in the flow field will scatter laser light that is directed upon them. A very sparse particle concentration, such as with

methanol, results in a very low data rate and small amplitude doppler bursts. A doppler signal of amplitude less than 200 millivolts will not be validated in the counter processors, resulting in a near zero data rate. The unusually long focal length of 16 feet also contributes to a smaller doppler signal in comparison to a short focal length, as the light scattered from the particles in the flow field travels a longer distance before being collected in the photomultiplier tube. The longer distance travelled by the scattered light results in reduced signal strength.

In this run, the doppler signal was boosted by adding small amounts of n-heptane to the fuel pan. This resulted in a more soot-laden flame of orange color. The flame was placed inside a wooden box with several small cut-out holes for ventilation. The box served to stabilize the flame from the air currents created by the nearby hood and air ducts, as well as to elevate the smoke particle concentration near the flame to boost the doppler signal strength.

Appendix C contains a listing of the data which has been grouped by doppler burst. Column #2 lists velocity in the vertical direction (meters per sec), column #3 lists velocity in the horizontal direction, and column #4 gives the sample interval time for that doppler burst from which the data is generated. For the vertical component, a positive number indicates the upward direction, while negative indicates downward. A positive number for the horizontal component indicates a direction to the right. After processing all of the data, the unweighted mean velocity, as calculated from all the data points, is shown along with the angle of the velocity vector, in reference to the vertical position.

V. CONCLUSIONS

In a continuing effort to understand the physical mixing processes in diffusion flames, NRL is currently involved in the development of a technique to measure the velocity profile in combustion processes using the laser doppler velocimeter. To date, the software has been developed for a computer-controlled LDV. Further work is required in the areas of seeding the flow field to boost the signal to noise ratio of the doppler signal to make this a viable research tool.

VI. REFERENCES

 Asalor, J. O. and Whitelaw, J. H., "The Influence of Combustion-Induced Particle Concentration Variations in Laser-Doppler Anemometry," Proceedings of the LDA Symposium, Copenhagen (1975).

- Stevenson, W. H., Thompson, H. D. and Gould, R. D., "Laser Velocimeter Measurements and Analysis in Turbulent Flows with Combustion, Part II," School of Mechanical Engineering, Purdue University, Report No. AFWAL-TR-82-2076, July 1983.
- Alexander, J. I., St. Aubin, H. J., Stone, J. P., Street, T. T. and Williams, F. W., "Large-Scale Pressurizable Fire Test Facility - FIRE I," NRL Formal Report 8643, 30 December 1982.

APPENDIX A - SOFTWARE DRIVER

```
Tabryr T=00003 IS ON CROO041 USING 00015 BLKS R=0000
 0001
       454B, R, L
             NAM DRYR,7
                           USER WRITTEN DRIVER(SELECT CODE=108)
 0002
0003
             EHT DRYR
 0004
             EXT ENTRISLIBRISLIBX
 0005
       * FIND ADDRESSES OF PASSED VARIABLES
 0006
 0007
       ICODE NOP
 ogo 8
       IDATA NOP
 0009
       IOL
             HOP
       ISCU
 0010
            NOP
 0011
       ISL
             NOP
 0012
       ICHT
             NOP
 0013
       IFLAG NOP
 00:4
       DRYR HOP
 0015
             JSB .ENTR
 0016
             DEF ICODE
 2017
 0018
       * DEFINE SELECT CODE OF 12566 BOARD
            EQU 108
       LDV
 0013
 0329
       . DETERMINE TASK:
 0 0 2 1
 0022
             LDA ICODE, I
                                 ICODE=1 DO INITIATION SECTION ONLY
             CPA THO
 0023
                                 ICODE=2 DO CONTINUATION SECTION ONLY
             JMP CONT
 0024
                                 ICODE=3 DO BOTH SECTIONS
 0025
 0026
       . INITIATION SECTION:
 0027
       **************
             LDB DC
                            KEEP 60000B IN B
 0028
 0029
             LDA KI
                            SAVE MEMORY SIZE TO CLEAR
0030
             CNA . INA
 0031
             STA CHT1
0032
             LDA ISL, I
                            USE ISCW LENGTH AS A COUNTER
 0033
             CHA, INA
             STA CHT2
LDA ISCU
 0034
0035
                            SAVE ADDRESS OF CONTROL BUFFER
0036
             STA PTRI
0037
 0038
       . DISABLE THE MEMORY PROTECT FENCE
0039
             JSB $LIBR
                            GO PRIVILEGED
 0040
             NOP
0041
       . CLEAR THE USED BUFFER MEMORY
0042
       CLEAR NOP
0043
 0044
             OTB LDV
 0045
             STC LDV.C
 0046
             CLC LDV
 0047
             ISZ CHT1
 0048
             JMP CLEAR
0049
       * SEND OUT CONTROL WORDS TO INTERFACE
0050
 0051
       SEND LDA PTRI, I
0052
             OTA LDY
 0053
             STC LDV, C
             CLC LDV
 0054
 0055
             ISZ PTRI
                          INCREMENT BUFFER ADDRESSES
 0056
             ISZ CHT2
                          AND COUNTER
 0057
             JMP SEND
 0053
```

```
0059
      . RETURN IF ICODE=1
            LDA ICODE,I
0060
0061
             CPA ONE
             JHP BYE
0062
0063
      . CONTINUATION SECTION :
0064
0065
0066
      CONT NOP
0067
             LDA ONE
                              ASSUME A SUCCESSFUL CALL
0068
             STA IFLAG, I
0069
            LDA IDL, I
                              USE BUFFER SIZE AS A COUNTER
0070
            CHA. INA
0071
             STA CHT3
0072
            LDA IDATA
                              STORE ADDRESS OF DATA BUFFER
0073
             STA PTR2
0074
                              KEEP 60000B IN B-REG
             LDB DC
             LDA ZERO
0075
                              INITIATE DATA COUNTER
0076
             STA ICHT, I
0077
0078
      . IF ICODE-3 NO NEED TO TURN OFF OPERATING SYSTEM
            LDA ICODE,I
CPA THRE
0079
0080
0081
             JMP MORE
0082
0083
      . TURN OFF OPERATING SYSTEM
0084
             JSB $LIBR
0085
             HOP
0086
0087
      * COLLECT DATA TRANSFER WORDS
0088
      HORE NOP
0089
             LDA THOU
0090
             STA WAIT
0091
             OTB LDV
             STC LDV.C
0092
0093
      LIST
            NOP
0094
             SFS LDV
0095
             JHP IHP
0096
             LIA LDV
0097
             CPA ZERO
                                PREVENT ZEROES FROM ENTERING
0098
             JHP HORE
                                BUFFER MEMORY
0099
             STA PTR2, I
0100
             ISZ PTR2
0101
             NOP
0102
             ISZ ICHT, I
0103
             NOP
0104
             ISZ CHT3
0105
             JHP HORE
0106
      . ENABLE THE MEMORY PROTECT FENCE & RETURN
0107
0108
      BYE
            NOP
0109
             CLC LDV
0110
             CLF LDV
0111
             JSB $LIBX
                              GO UNPRIVILEGED
0112
             DEF ++1
0113
             DEF ++1
             JMP DRYR, I
0114
0115
      * ROUTINE TO WAIT THOU TIMES FOR FLAG
0116
0117
      IMP
            HOP
9118
             ISZ WAIT
```

Ł

```
JMP LIST
0119
0120
0121
             STA IFLAG, I
             CLC LDV
CLF LDV
JSB $LIBX
0122
0123
0124
             DEF ++1
DEF ++1
JMP DRYR, I
0125
0126
0127
0128
0129
0130 . VARIABLE STORAGE
0131
      DC
             OCT 60000
0132
      K 1
             DEC 512
0133
      ZERO
             DEC 0
0134
      OHE
             DEC 1
             DEC 2
0135
      TUO
0136
      THRE
             DEC 3
      THOU
             DEC -32000
0137
      CHTI
             NOP
0138
0139
      CHT2
             HOP
0140
      CHT3
             HOP
             HOP
0141
      PTRI
0142
      PTR2
             HOP
             NOP
0143
      VAIT
0144
              EHD
```

APPENDIX B - MAIN FORTRAN PROGRAM

```
. SLDV
       T=00003 IS ON CR00041 USING 00094 BLKS R=0000
      FTN4,L
0001
0002
0003
      C
0004
            PROGRAM LDY(3,91), REV 830720- TO USE THE LDY
0005
0006
            WRITTEN BY STEVE LUSTIG AND CAROLYN KAPLAN
0007
0002
      C----FUNCTIONS:------
0009
      C 1) TO TEST THE COINCIDENCE FILTER BOARD
      C 2) TO OPERATE THE DISA LASER DOPPLER AMENOMETER
0010
           SYSTEM HARDWARE FROM HP-21MX SOFTWARE
0011
0012
        3) TO COLLECT DATA FROM DISA LASER DOPPLER ANEMOMETER
0013
           BUFFER INTERFACE
      C 4) TO CREATE DISK/TAPE FILES OF THIS DATA
0014
0015
      C 5> TO ASSEMBRE DATA & WRITE OUT ON LP (LU6)
      C 6) TO COMPUTE REYNOLD'S STRESS CORRELATIONS (& ETC)
0016
           FROM ASSEMBLED DATA
0017
0012
      C---- MISCELLANEOUS-----
0019
      C 1> THE DRIVER OPERATING THE LDV IS CURRENTLY A
0020
0021
           RELOCATABLE VERSION WHICH IS LOADED WITH THIS PROGRAM.
0022
           THE VARIABLES PASSED TO THE DRIVER(LDRVR) MUST NOT BE
0023
           CHANGED .
0024
      C
        2) THE LAST WORD IN ARRAY SCW MUST BE THE 'STATIC CONTROL C
           WORD' WHICH IS USED TO 'TURN ON' DESIRED CHANNELS. PRE C
0025
           CEEDING CONTROL WORDS OPERATE PERIPHERAL BOARDS/DEV.S.
0026
0027
      C
        3> THE SORTING ROUTINE [MODE=3] WILL MOST LIKELY HAVE TO
0028
           BE MODIFIED AS NORE HARDWARE IS ADDED TO THE LOV SYS.
          RUNNING THIS PROGRAM WITH THE FTN4.0 OPTION WILL ACTI- C
0029
0030
           VATE PRINTING SUBROUTINES IN TASK 3.
      C
0031
      0032
0033
            COMMON LU(5)
0034
            INTEGER ISCU(10), IDATA(512), IFLAG, TITLE(20), IDCB(656)
            INTEGER TAPELU, ITINE(5), IBUFF(20), IDCB2(144), NAME(3)
0035
0036
            INTEGER NAME2(3), PLOC, PARTCH, STIME, ICHT, WAIT
            REAL FTI. FTW. FTC. TTI. DT. COUNT. FORM1. FORM2. RD ATA
0037
0038
            DATA TITLE/40+2H
0039
      C---- DETERMINE THE LU OF THE DEVICE------
0040
0041
            CALL RMPAR(LU)
0042
0043
0044
      C----DETERMINE TASK-------
0045
0046
            WRITE(LU,3)
        2
0047
            FORMAT(1X, 'ENTER:
                                1
                                    TO TEST LDV BY TAKING SAMPLES'
0048
           +/17X,'(WRITE TO ANY NON-DISC LU)'
           +/11X,'2
0049
                      TO COLLECT DATA FROM THE LDV INTERFACE'
0050
           +/17X,'(COLLECT FOR A GIVEN AMOUNT OF TIME,
           +/17X, WRITE RAW DATA TO A CREATED DISC FILE)
0051
                      TO SORT AND ASSEMBLE COLLECTED DATA
0052
           +/17X,'(READ RAW DATA FROM ANY DISC LU,'
0053
0054
           +/17X, 'ASSEMBLE DATA, WRITE TO DISC FILE)'
                      TO COMPUTE TURBULENCE STRESS CORRELATION'
0055
           +/11X,'4
           +/17X, 'VALUES FROM ASSEMBLED DATA FILES'
0056
                      TO END THIS PROGRAM'//11X,'_')
0057
           +/11X,15
0058
            READ(LU,+) INODE
```

```
G0 T0 (100,200,300,400,9999), IMODE
0059
0060
       100 CONTINUE
0061
      C----TEST THE SYSTEM AND WRITE TO ANY NON-DISC LU------C
0062
0063
0064
      C ASK FOR AMOUNT OF DATA AND OUTPUT LU
0065
             ISL = 0
0066
             IDL=512
0067
             ICOUNT=0
            WRITE(LU, 109)
0068
            FORMAT(' ENTER # DATA WORDS TO COLLECT((=512)(151: _')
0069
            READ(LU, +)IDL
0070
0071
            IF( IDL GT . 512 ) IDL=512
9072
            URITE(LU, 130)
            FORMAT(' ENTER OUTPUT LU(12): _')
0073
       130
            READ(LU, 131)IOUT
0074
0075
            FORMAT(12)
       131
0076
0077
      C READ STATIC CONTROL WORDS AND ZERO BUFFER
0078
            WRITE(LU,103)
0079
       103
           FORMAT(' ENTER A CONTROL WORD (000000 TO STOP) [06]: _')
0080
            READ(LU, 104) IEHTER
0081
       104
            FORMAT(06)
0082
            IF(IENTER.EQ. 0B) GO TO 105
0083
            ISL = ISL+1
0084
            ISCW(ISL)=IENTER
0085
            GO TO 102
0086
       1.05
            DO 110 I=1.512
0087
            IDATA(I)=0
       110
0088
            IF(ISL.EQ.0) G0 T0 2
0089
      C
0090
            ICODE=3
            CALL DRYR(ICODE, IDATA, ISCN, ISL, ICHT, IFLAG)
0091
0092
      C
      S WRITE OUT ANY DATA COLLECTED
0093
0094
            IF(ICHT.GT.O) WRITE(IOUT,108) (IDATA(I),I=1,ICHT)
0095
       108
            FORMAT(65(/8(1X,06)))
0096
            WRITE(IOUT, 120) ICHT, IFLAG
            FORMAT(' ICHT = ', 13/' IFLAG = ', 12//)
0097
       120
0098
0099
      C COLLECTED ENOUGH ?
            ICOUNT=ICOUNT+ICHT
0100
            IF(ICOUNT.LT.IDL) GO TO 105
0101
0102
            GO TO 2
0103
      C
0104
       200 CONTINUE
     C----COLLECT DATA FROM LDV INTERFACE-----
0105
0106
      C SELECTIONS OF STATIC CONTROL WORDS CAN BE ALTERED/STORED HERE C
      C FOR ROUTINE DATA COLLECTIONS. THE DISA MANUALS SHOULD BE
                                                                         C
0107
0108
      C CONSULTED BEFORE MAKING HARDWARE CHANGES ANY HARDWARE MODI-
                                                                         C
      C FICATION WILL PROBABLY RESULT IN THE NEED TO HODIFY CONTROL
0109
                                                                         C
     C WORDS. THIS IS A GOOD PLACE TO STORE CONTROL WORDS!
                                                                         C.
0110
0111
            ISCU(1)=010011B
0112
            ISCW(2)=0223428
            ISCU(3)=0000548
0113
0114
            1 SL = 3
            IDL = 512
0115
0116
            MAIT=5
0117
     C 1. ISCW(1) TURNS OH THE COINCIDENCE FILTER BOARD AND CHANNELS C
0118
           ONE AND TWO.
```

```
0119
      C 2. ISCU(2) SETS A 125 MICROSEC TIME INTERVAL FOR A 0.1
           HICROSEC TIME BASE FOR THE COINCIDENCE FILTER BOARD.
0120
           ISCU(3) SAMPLES BOTH PART CHANNELS OF CHANNELS 1,2 AND
0121
                                                                           £.
            CORRESPONDING TO PL, P2 AND P4 ON THE INPUT MULTIPLEXER
0122
                                                                           C
0123
            BOARD
       4. THE DRIVER ASSUMES THAT THE ADDRESS TO DROP REQUEST B IS
0124
      C
                                                                           C
0125
           SET AS 68 (THE SWITCH ON THE 57GXXX BORD IS ON 6)
0126
        TAKE DATA CALLING USER WRITTEN DRIVER AND PASS THE FOLLOWING
                                                                           C
0127
      C
        VARIABLES:
0128
      C ICODE
                      : =1 PERFORM INITIATION SECTION OF DRIVER ONLY
                                                                           С
                        (CLEAR BUFFER MEMORY AND SEND CONTROL WORDS)
0129
                                                                           C
0130
                        =2 PERFORM CONTINUATION SECTION ONLY
0131
                        (COLLECT ANY DATA IN BUFFER MEMORY)
                        =3 DO BOTH SECTIONS WITHOUT RETURNING TO THIS
0132
0133
                            PROGRAM IN BETWEEN
0134
        IDATA(512)
                      : INTEGER ARRAY OF DATA TRANSFER WORDS COLLECTED
0135
       IDL
                      : INTEGER HUMBER OF DATA SAMPLES TO COLLECT( = <512)
      C
0136
      C
       ISCU(10)
                      : INTEGER BUFFER OF CONTROL WORDS TO SEND TO THE C
0137
                         57G20 CONTROL SYSTEM. THE LAST WORD SHOULD BE
                                                                           C
                        THE STATIC CONTROL WORD (SCW)
0138
                      : INTEGER HUMBER OF CONTROL WORDS TO BE SENT.
0139
      C ISL
                                                                           C
                      : # OF DATA WORDS COLLECTED BY THE DRIVER
0140
      C ICHT
                     (INITIATED AT THE BEGINNING OF EVERY CALL) = 1 MEANS THAT THE DRIVER COLLECTED DATA NORMALLY.
0141
0142
        IFLAG
                     =2 THE DRIVER HAD TO WAIT TOO LONG FOR THE FLAG
0143
0144
                        FLIP-FLOP TO SET WHILE RTE WAS SUSPENDED
0145
      C
        NOTE THAT IDATA BUFFER IS NOT CODED ON DISC IN ASCII CHARCTR
      C SET. THUS IT IS IMPOSSIBLE TO SEE THIS DATA W/OUT SOFTWARE.
0146
0147
      C IN ADDITION, THE DATA BUFFER IS NOT ZEROED BETWEEN DRVR CALLS C
        TO SAVE TIME.
                       THIS PROGRAM KNOWS HOW MANY NEW WORDS HAVE BEEN C
0148
0149
        COLLECTED BY THE COUNTER ICHT. THUS ONLY THE FIRST ICHT WORDS C
0150
      C IN BUFFER IDATA HAVE BEEN COLLECTED AFTER ANY DRVR CALL!
0151
      C---
0152
            COUNT=0.0
            DATA TITLE/40+2H
0153
0154
      C GET INFORMATION FROM USER
0155
0156
             URITE(LU, 211)
            FORMAT(5X,' ENTER CARTRIDGE #, OUTPUT FILE NAME(12, A6) _')
0157
0158
            READ(LU, 202) ICR, NAME
0159
            FORMAT(12,1X,3A2)
       202
0160
            WRITE(LU, 204)
0161
            FORMAT(5x,' ENTER APPROXIMATE TIME IN'
            +/6X.' SECONDS TO COLLECT DATA([10]: _')
0162
0163
            READ(LU, 208) STIME
0164
       208
            FORMAT(I10)
0165
             WRITE(LU, 209)
            FORMAT(5X,' ENTER A TITLE[A36]')
0166
       209
0167
             READ(LU, 210)(TITLE(I), I=2, 20)
0168
       210
            FORMAT(2042)
0169
      C
0170
      C CREATE AND OPEN A TYPE 3 SEQ. ACCESS FILE
0171
             ISIZE=-1
0172
             IOPTH=0
0173
             ITYP=3
0174
             ISECU=0
0175
             IDCBS=640
0176
            CALL CREAT(IDCB, IERR, NAME, ISIZE, ITYP, ISECU, ICR, IDCBS)
0177
             IF(IERR.LT.0) CALL ERR(1, IEPR,1)
             CALL OPEN(IDCB, IERR, NAME, IOPTH, ISECU, ICR, IDCBS)
0178
```

)

)

```
0179
             IF(IERR LT.0) CALL ERR(2, IERR, 2)
0180
      C
0181
      C WRITE OUT A TITLE
0182
            NWORDS=20
             CALL WRITF(IDCB, IERR, TITLE, NWORDS)
0183
0184
             IF(IERR.LT.0) CALL ERR(5,IERR,3)
0185
      C GET STARTING TIME
0186
0187
            CALL EXEC(11.ITIME, IYEAR)
0188
            FTI=(ITIME(4)+60+ITIME(3))+60+ITIME(2)
0189
0190
      C COLLECT AND WRITE OUT A BUFFER OF DATA
0191
        203 ICODE=1
0192
             CALL DRVR(ICODE, IDATA, IDL, ISCW, ISL, ICNT, IFLAG)
0193
                  CALL EXEC(11, ITIME, IYEAR)
0194
                  TTU=(ITIME(4)+60+ITIME(3))+60+ITIME(2)
0195
        205
                       CALL EXEC(11, ITIME, IYEAR)
0196
                        NTU=WAIT+(ITIME(4)+60+ITIME(3))+ITIME(2)
0197
                        IF(HTW.LT.TTW) GO TO 205
0198
             ICODE=2
0199
             CALL DRYR(ICODE, IDATA, IDL, ISCN, ISL, ICNT, IFLAG)
0200
0201
         WRITE BUFFER TO TERMINAL, IF DESIRED
0202
      D
             IF(ICHT.GT.0) WRITE(LU,108)(IDATA(I),I=1,ICHT)
             WRITE(LU, 120) ICHT, IFLAG
0203
      D
0204
      C
0205
      C
         WRITE OUT DATA TO DISC FILE
             IF(ICHT.GT.O) CALL WRITF(IDCB, IERR, IDATA, ICHT)
IF(IERR.LT.O) CALL ERR(5, IERR, 5)
0206
0207
             COUNT=COUNT+ICHT
0208
      C
0209
      C CHECK TIME
0210
            CALL EXEC(11.ITIME.IYEAR)
0211
             TTI=(ITIME(4)+60+ITIME(3))+60+ITIME(2)
0212
            DT=TTI-FT1
0213
             IF(DT.LT.STIME) GO TO 203
0214
0215
             RATE = COUNT/DT
0216
0217
      C WRITE OUT . SAMPLES AND CLOSE THE DISC FILE
             WRITE(LU, 207) COUNT, RATE, NAME, ICR
0218
           FORMAT(10X, 'SAMPLING TIME OVER: '
0219
            +/10X,G15.3, DATA WORDS COLLECTED //7X, AT /,G15.3,
0220
            +' WORDS/SEC'/7X,'OH ',3A2,'::',12//)
0221
             CALL LOCF(IDCB, IERR, IREC, IRB, IOFF, JSEC)
0222
0223
             IF(IERR.LT.O) CALL ERR(3, IERR,6)
0224
             ITRUN=JSEC/2-IRB-1
0225
             CALL CLOSE(IDCB, IERR, ITRUH)
0226
             IF(IERR.LT.O) CALL ERR(4,IERR,7)
0227
      C GO BACK TO COMMAND SECTION
0228
0229
             GO TO 2
0230
      C-----SORT & REDUCE RAW DATA FILES-----
0231
       300 CONTINUE
0232
0233
0234
      C PREPARE INPUT DISC FILE WITH RAW DATA
0235
             IOPTH=0
0236
             ITYP=3
             ISIZE=-1
0237
             I SE CU= 0
0238
```

```
1239
             IDCBS=640
0240
             IDC852=128
             WRITE(LU, 304)
0241
0242
            FORMAT(5x, 'ENTER CARTRIDGE # , INPUT FILE NAME [12, 46] _ ')
0243
             READ(LU, 202) ICR, NAME
0244
             CALL OPEN(IDCB, IERR, NAME, IOPTH, ISECU, ICR, IDCBS)
0245
             IF(IERR.LT.0) CALL ERR(2,IERR,8)
0246
0247
      C CREATE AND OPEN OUTPUT FILE
0248
            WRITE(LU, 211)
0249
             READ(LU, 202) IOUT, NAME 2
0250
            CALL CREAT(IDCB2, IERR, NAME2, ISIZE, ITYP, ISECU, IDUT, IDCBS2)
0251
             IF(IERR, LT. 0) CALL ERR(1, IERR, 9)
0252
             CALL OPEN(IDCB2, IERR, NAME2, IOPTH, ISECU, IOUT, IDCBS2)
0253
             IF(IERR.LT.0) CALL ERR(2, IERR, 10)
0254
0255
      C WRITE ASSEMBLED DATA TO A NON-DISC LU, IF DESIRED
0256
      D
            WRITE(LU, 305)
0257
      0305
            FORMAT(5x, ENTER A NON-DISC LU TO SEE ASSEMBLED DATA:
0258
             READ(LU, 131) ILU
0259
             WRITE(LU, 306) ILU
0260
      0306
            FORMAT(5X, DO YOU HEED TO CONTROL LU ',[2,' ?[Y/N] _')
             READ(LU, 307) IANS
0261
      0307
0262
            FORMAT(A1)
0263
      D
             IF(IANS.EQ. 1HY) CALL TAPE(ILU)
0264
0265
      C READ IN TITLE AND WRITE TO OUTPUT FILE
            HUORDS=20
0266
0267
             CALL READF(IDCB, IERR, IBUFF, NWORDS, ILEN)
0268
             IF(IERR.LT.O) CALL ERR(4,IERR,11)
             CALL WRITF(IDCB2, IERR, IBUFF, NWORDS)
0269
0270
             IF(IERR.LT.0) CALL ERR(5, IERR, 12)
0271
            HCHRS=40
0272
             CALL CODE(NCHRS)
             READ(IBUFF, 210)(TITLE(I), I=1, 20)
0273
0274
0275
      C INITIATE STATISTICAL QUANTITIES, IPASS=1
0276
            IPASS=1
0277
             CALL STAT(IPASS, ICHTR, RDATA, ITYPE)
0278
      C
0279
      C CALL FIRST PRINTING ROUTINE, IF DESIRED
0280
      D
            CALL PRINICTITLE, ILU)
0281
0282
      C READ IN A RECORD OF DATA OF UP TO 512 WORDS
0283
            HCHRS=22
0284
             HWRDS=11
0285
             HWORDS=512
0286
             IPASS=2
0287
        309 CONTINUE
0288
            DO 308 I=1,512
0289
        308 IDATA(1)=0
0290
             CALL READF(IDCB, IERR, IDATA, NUORDS, ILEN)
0291
             IF(IERR.LT.0) CALL ERR(6, IERR, 14)
0292
             IF(ILEN.LT.0) GO TO 311
0293
0294
      C PROCESS THE DATA ON A RECORD
0295
            ISKIP=0
0296
            DO 310 I=1, ILEH
0297
      C
0298
      C
                 ASSEMBLE THE DTW INFORMATION
```

2

```
CALL ASM8L(IDATA(I), PARTCH, PLOC, FORM1, FORM2, IBO, ITO)
0299
0300
0301
                SORT THE INFORMATION IN THE DTW
              CALL SORT(IDATA(I), PARTCH, PLOC, ITYPE, ICHTR, FORM1, FORM2, RDATA)
0302
0303
              IF(ITYPE.EQ 6) GO TO 310
0304
                PRINT/WRITE THE DATA, IF DESIRED
0305
0306
              CALL PRINZ(ICHTR, ITYPE, RDATA, IBO, ITO, ILU)
0307
0308
                DO STATISTICAL SUMMATIONS, IPASS=2
0309
              CALL STAT(IPASS, ICHTR, RDATA, ITYPE)
9319
0311
                GROUP VELOCITY & INTRVL DATA (PER BURST)& WRITE TO OUTPUT FILE
0312
              CALL GATHR(ISKIP, ICHTR, ITYPE, RDATA, IDCB2)
9313
       310 CONTINUE
2314
     C READ THE NEXT RECORD IN THE FILE
0315
9316
            GO TO 309
0317
      £.
0318
      C DO FINAL STATISTICAL COMPUTATIONS ON A BUFFER OF DATA, IPASS=3
       311 IPASS=3
0319
0320
            CALL STAT(IPASS, ICHTR, RDATA, ITYPE)
0321
0322
      C CLOSE THE BUTPUT FILE
0323
            CALL LOCF(IDCB2, IERR, IREC, IRB, IOFF, JSEC)
0324
            IF(IERR.LT.0) CALL ERR(3,IERR,16)
0325
            ITRUH=JSEC/2-IRB-1
            CALL CLOSE(IDCB2, IERR, ITRUN)
0326
            IF(IERR.LT.0) CALL ERR(4, IERR, 17)
0327
0328
0329
     C PLACE EOF ON ILU IF IT IS A TAPE
0330
            ICHUD=ILU+100B
            CALL EXEC(3, ICHND, 0)
0331
0332
     C RETURN TO COMMAND SECTION
0333
0334
           G0 T0 2
0335
0336
       400 CONTINUE
0337
      C----COMPUTE REYNOLDS STRESSES-----C
       9999 WRITE(LU, 9998)
0338
0339
       9998 FORMAT(10X, '**END OF LDV**')
0340
            END
0341
            SUBROUTINE SORT( IDATA, PARTCH, PLOC, ITYPE, ICHTR, FORM1,
0342
           +FORM2, RDATA), TO DECIPHER DTW'S
0343
0344
            COMMON LU(5)
0345
            INTEGER P1234(4), DIGIO(4), MODE(4), PARTCH, PLOC
            REAL FSHIFT(4), LAMBDA(4), THETA(4), FORM1, FORM2, TBASE, RDATA
0346
0347
     C**** HOTICE TO USER MAKING MODIFICATIONS *****
0348
0349
      C INDICATE DEVICES CONNECTED TO THE 57G120 INPUT MULTIPLEXER:
0350
        THE USER TELLS THIS PROGRAM HOW THE 57G20 CABLES ARE CON-
0351
           HECTED BY SETTING THE MEMBERS OF THE FOLLOWING ARRAYS:
0352
           P1234(X)=Y WHERE
                                X=1,2,3, OR 4 CORRESPONDING TO PINS
                                  P1, P2, P3, P4 RESPECTIVELY.
0353
                                                                         C
                                Y=1,2 FOR COUNTER WHICH WILL BE REF-
0354
                                      EREHCED AS 1.2
0355
0356
                                Y=9 FOR THE COINCIDENCE FILTER BOARD
0357
                                Y=0 IF NOTHING IS CONNECTED AT ALL
9358
                                  OR FOR THE COINCIDENCE FILTER BOARD
```

```
0359
     C
           DIGIO(X)=A WHERE
                               X IS THE SAME AS ABOVE
0360
                               A=1 IF THE DIGITAL I/O SLOT #1 CORRES
0361
                                 TO PIN X OR A=2 IF DIGITAL I/O SLOT
0362
     C
                                 #2 APPLIES
                               A=O IF A COUNTER IS NOT CONNECTED
0363
      C
0364
           FSHIFT(X)=B WHERE
                               X IS THE SAME AS ABOVE
                                                                       C
                               8 IS THE FREQUENCY SHIFT AS SHOWN ON
0365
                                                                       C
0366
                                 THE 55N14 FREQ GENERATOR, I E B IS
                                                                       C
0362
                                 POSITIVE IF THE >40MHZ SIDE IS USED
                                                                       C
                                 AND NEGATIVE IF (40MHZ B IS ONLY THEC
0368
0369
      C
                                 VALUE SET ON THE 55H14 PANEL
0370
      C
           MODE(X)=C
                        WHERE
                               X IS THE SAME AS ABOVE
0371
                               C=1 IF DATA IS COLLECTED IN THE FIXED
                                                                       C
0372
                                 MODE ON THE COUNTER PROCESSOR,
                                                                       C
0373
                               C=2 IF THE COMBINED HODE IS USED.
                               C=3 IF THE VARIABLE FRINGE MODE, AND
0374
      C
0375
                               C=4 IF THE TRANSIT MODE IS USED
                               C=9 IF HO COUNTER IS ON X.
      С
0376
                                                                       C
0377
      C
            TBASE
                              IS THE FLOATING POINT CLOCK TIME BASE
                              OH JUMPER S2 OF THE CFB. IT MAY BE 0.1 C
0378
      C
0379
                              1, OR 10 USEC
                         WHERE X IS THE SAME AS ABOVE
0380
      £
            LAMDA(X)=D
0381
      C
                               D IS THE WAVELENGTH OF LASER EHR USED
                               ON CHANNEL X.
0382
      C
0383
            THETA(X)=E WHERE
                               X IS THE SAME AS ABOVE.
0384
     C
                               E IS THE ANGLE BETW/ THE CORRES. BEAMS C
0385
      C
                                 ASSOC W/ THIS COUNTER (FOR COAXIAL)
                                                                       C
                                 BEAMS THIS IS ALWAYS CONSTANT)
0386
      C
0387
0388
            DATA P1234/2,1.0,9/
0389
            DATA DIGIO/2,2,0,0/
0390
            DATA FSHIFT/-10 0E+03,-10.0E+03,0.0,0.0/
0391
            DATA MODE/1,1,9,9/
0392
            DATA TBASE/0.1E-06/
            DATA LAMBDA/ 488.0E-09.514.5E-09.0.0.0
0393
            DATA THETA/0.019330253,0.019701965,0.0,0.0/
0394
0395
0396
     C------
0397
     C ITYPE IDENTIFIES THE TYPE OF INFORMATION THAT A DATUM GIVES: C
     C ITYPE= 1 FOR DOPPLER FREQUENCY DATA
0398
2399
               2 FOR SAMPLE INTERVAL DATA,
0400
               3 FOR NUMBER OF FRINGES,
               1 FUR BURST TIME DATA,
0401
      C
               5 FOR VELOCITY DATA, AND
0402
                                                                      C
0403
              6 IDATA=0
                                                                      C
0494
0405
     C IDENTIFY THE TYPE OF DATA
0406
            ITYPE=DIGIO(PLOC)*PARTCH
0407
            IFCITYPE.EQ.2 AND.PARTCH.EQ.18) ITYPE=3
0408
0409
      C SAMPLE INTERVAL HAS ONLY ONE PART CHANNEL OF INFO:
3410
            IF(P1234(PLGC) EQ.9 AND.PARTCH.EQ 18) ITYPE=2
0411
            IF(P1234(PLOC) EG 9 AND PARTON EG 28) ITYPE=6
0412
            IFCIDATA . EQ. O > ITYPE=6
0413
0414
      C CONVERT DIGITAL OUTPUT TO DOPPLER FREQ, BURST/TRANSIT TIME,
0415
      C NUMBER OF FRINGES, OR SAMPLE INTERVAL TIME. THIS IS DONE AS
     C PRESCRIBED ON P27 OF DISA COUNTER INSTRUCTION MANUAL. ALSO
0416
     C SUBTRACT FREQUENCY SHIFTING AND CONVERT TO FREQUENCY TO MELOCITY
0417
            IF(MODE(PLOC).GT.2 + GO TO 353
0418
```

```
IFCITYPE EO () RONTH=FORM1+14 96+FSHIFT(PLOC)
0413
               IF(ITYPE EQ 4) ROATA=16.E+09/FORM1+FSHIFT(PLOC)
11420
9421
               IF(ITYPE.EQ 4) ITYPE=1
942.
            IF MODE(PLOC) LE 2) GO TO 354
                IFCITYPE.EQ.1) RDATA=0.5348/FORH1+FSHIFT(PLOC)
0423
                IF(1TYPE EQ.4) RDATA=0.5E-09*FORM1+FSHIFT(PLOC)
0424
9425
               IF(ITYPE EQ 4) ITYPE=1
0426
      C CONVERT DOPPLER FREQUENCY TO VELOCITY
0427
       354
           IF(1TYPE NE 1) GO TO 355
            IF(1.EQ.1) GO TO 355
0428
            FACTOR=LAMBDA(PLOC)/(2.*SIN(THETA(PLOC)/2.))
0429
9439
            RUNTH= RUNTH+FACTOR
0431
            ITYPE=5
       355
           IF(ITYPE.EQ.2) RDATA= FORM1*TBASE
0432
            IF(ITYPE EQ 3) RDATA= FORM2
0433
0434
      C IDENTIFY WHICH COUNTER
0435
            ICNTR=P1234(PLOC)
U436
9437
0438
      C LEHVE
0439
            RETURN
0440
            EHD
0441
            SUBROUTINE GATHR(ISKIP, ICHTR · ITYPE, RDATA, IDCB2), TO GROUP DATA BY
0442
0443
           + BURSTS
0444
            COMMON LU(5)
0445
            DIMENSION IBUFF(33), ITYP(10), ICHT(10), RDAT(10)
0446
            IADDRS=IENTR*ITYPE
9447
0448
         IF THIS IS THE FIRST CALL, SAVE THE FIRST ADDRESS
            IF( ISKIP, EQ. 0) G0 T0 506
0449
0450
0451
         IF IT IS NOT THE FIRST CALL, CHECK TO SEE IF ADDRESS EQUALS
      C
         FIRST ADDRESS
0452
      Ç
9453
            IF(IADDRS.NE.IAD)GO TO 522
0454
0455
         IF ADDRESS EQUALS FIRST ADDRESS, ZERO THE ARRAY, AND RESET IPTR
9456
        521 IPTR=0
9457
            DO 508 [=1,3
0458
            ITYP(I)=0
0459
            ICHT(I)=0
            RDAT(I)=0.0
2460
0461
        508 CONTINUE
0462
      C
          IF DATA IS VELOCITY OR SAMPLE INTERVAL DATA WRITE IN ARRAY
0463
      С
        522 IF(ITYPE HE.2 AND. ITYPE. NE 5) RETURN
0464
            IPTR=IPTR+1
0465
0466
            ISKIP#1
0467
            IDUM = ICHTR
J468
            IF(ITYPE.EQ.2)IDUM=3
0469
            ITYP(IDUM)=ITYPE
0470
            ICHT(IDUM)=ICHTR
0471
            RDAT(IDUM)=RDATA
9472
     C
         IF ARRAY IS FULL. WRITE ARRAY TO BUFFER
9473
0474
            IFCIPTR.NE.35 RETURN
0475
            CALL CODE(NORRS)
0476
            WRITE(IBUFF,501)(ICHT(IPTR),ITYP(IPTR),RDAT(IPTR),IPTR=1,3)
        501 FORMAT(3(1X,12,1X,12,1X,G15.8))
9477
0478
```

```
0479
      C
          WRITE DATA IN BUFFER TO OUTPUT FILE & RETURN
0480
            CALL WRITF(IDCB2, IERR, IBUFF, NWRDS)
            IF(IERR.LT.0)CALL ERR(5,IERR,15)
0481
0482
            RETURN
0483
      C
0484
         IF FIRST CALL, SAVE FIRST ADDRESS IN VARIABLE IAD, AND
0485
      C
         THEN ZERO THE BUFFER
0486
        506 IAD=IADRS
            NCHRS=66
0487
0488
            HWRDS=33
0489
            GO TO 521
0490
            END
0491
0492
            SUBROUTINE STAT(IPASS.ICHTR.DATA, ITYPE). TO DO STATISTICS
0493
            COMMON LU(5)
0494
            REAL MAG, ANGL, SUM1, SUM2, AVG1, AVG2
0495
            INTEGER N1, N2
0496
            GO TO (100,200,300), IPASS
0497
0498
      C INITIATE QUANTITIES
0499
       100 CONTINUE
0500
            N1=0
0501
            N2=0
0502
            SUM 1 = 0
0503
            SUM 2 = 0
0504
            RETURN
0505
      C
0506
      C DO SUMMATIONS, ETC
0507
       200
            CONTINUE
0508
             IF(ITYPE: HE.5) RETURN
0509
            IF(ICHTR.NE.1) GB TG 201
0510
            N1=H1+1
0511
            SUM1=SUM1+DATA
            RETURN
0512
0513
            H2=H2+1
            SUM2=SUM2+DATA
0514
0515
            RETURN
0516
0517
      C FINISH CALCULATIONS
       300 CONTINUE
0518
0519
            AVG1=SUN1/N1
            AVG2=SUM2/N2
0520
            MAG=SQRT(AYG1+AYG1+AYG2+AYG2)
9521
0522
            AHGL#ATAH(AVG2/AVG1)
0523
            ANGL = ANGL + 180 / 3. 141592654
            WRITE(6,301) AYG1, AYG2, MAG, AHGL
0524
0525
            WRITE(LU, 301) AVG1, AVG2, NAG, ANGL
0526
       301 FORMAT(//' UNWEIGHTED MEAN VELOCITIES:
            +/' <CHTR 1> ',G15.8/' <CHTR 2> ',G15.8
0527
           +//"
                 RESULTANT VECTOR: '/' MAGNITUDE = ',G15.8
0528
            +/' INCLUDED ANGLE = ',G15.8//)
0529
            RETURN
0530
0531
            END
0532
                                                                        ASSEMBLE
0533
            SUBROUTINE ASMBL(DTW, PARTCH, PLOC, FORM1, FORM2, 180, 170),
0534
            + DATA FROM DTW'S
      C THE FORMAT OF P.D. & T TYPE DATA:
0535
0536
         + DTC + BO CC1 CC2 + M7 E3 E2 + E1 E0 M6 + M5 M4 M3 + M2 M1 M0 +
0537
      C
0538
```

```
S THE FORMAT OF FRINGE DATA:
3539
0540
           0541
        + DTC + 80 CC1 CC2 +
                             F7 F6 + F5 F4 F3 + F2 F1 F0 +
0542
        9543
          COMMON LU(5)
0544
          INTEGER DIW, PARTCH, PLOC, IBO, ITO, IMAN, IEXP
V545
          REAL FORM1, FORM2
0546
          PARTCH=IAND(DTW, 100000B)/1000008+18
0547
          PLOC=IAND(DTW,030000B)/10000B+1B
0548
           IBO = IAND(DTW.40000B)/40000B
0549
          ITG=IAND: DTW, 40008 )/40008
          IMAH = [AND (DYG - 1778) + ITO + 2008
0550
          1EXP=14ND(DTW.36008)/2008
0551
0552
          FORM1=IMAH+(2. ++IEXP)
0553
          FORM2=IAND(DTW,37708)/108
          RETURN
0554
9555
          END
9556
     0557
          SUBROUTINE PRINI(TITLE, ILU), FIRST PRINTING ROUTINE
0558
          COMMON LU(5)
0559
          INTEGER ISCH(10), TITLE(20), STARS(40)
9569
          DATA STARS/2H .39*2H**/
9561
          WRITE(ILU,11> (TITLE(I), I=1,20), (STARS(I), I=1,40)
0562
          FORMAT(1H1, 20A2/1X, / COUNTER #/,
      11
          +6x.'PART CHANNEL', 15x, 'DATA SENT', 5x, 'OVFLOW BIT', 3x, 'TIME GUT'
0563
0564
          +/40A2)
0565
          RETURN
V566
     0567
          SUBROUTINE PRINZ(CNTR, ITYPE, RDATA, IBO, ITO, ILU), SECOND PRINT ROUT:
0568
0569
          + I HE
0570
          COMMON LU(5)
9571
          INTEGER CHTR, PARTCH, 180, ITO, ALPHA(15,5)
0572
          REAL RDATA
0573
          DATA ALPHA/2HDO, 2HPP, 2HLE, 2HR , 2HFR, 2HEQ, 2HUE, 2HNC, 2HY , 6+2H
0574
          +2HSA,2HMP,2HLE,2H I,2HNT,2HER,2HVA,2HL ,7+2H
0575
          +2H# ,2H0F,2H F,2HRI,2HNG,2HES,9+2H ,
9576
          +2HBU,2HRS,2HT ,2HTI,2HME,10+2H
0577
          +2HYE,2HL0,2HCI,2HTY,11*2H
     C DO NOT PRINT THE SAME DATUM OVER & OVER
0578
0579
          IF(RDATA.EQ.PDATA) RETURN
0580
     C
          PDATA=RDATA
0531
          WRITE(ILU, 10) CHTR, (ALPHA(I, ITYPE), I=1, 15), RDATA, IBO, ITO
0582
      1.0
          FORMAT(5%, I1, 8%, 15A2, G15, 4, 9%, I1, 8%, I1)
0533
          RETURN
9584
          END
0585
     0536
          SUBROUTIHE ERR(IPASS IERR, ILOC), TO FIND DISC & SOFTWARE ERRORS
0587
          COMMON LU(5)
0588
0589
     C
0590
           IF(IPASS.EQ.1) WRITE(LU,1)IERR, ILOC
0591
          IF(IPASS.EQ.2) WRITE(LU,2)IERR, ILOC
0592
           IF(IPASS.EQ.3) WRITE(LU,3) IERR, ILOC
0593
           IF(IPASS.EQ.4) WRITE(LU,4) IERR, ILOC
0594
           IF(IPASS.EQ.5) WRITE(LU,5)IERR, ILOC
0595
           IF(IPASS.EG.6) WRITE(LU,6) IERR, ILOC
9536
           IF((ERR GE.O) RETURN
0537
          3708
0538 C
```

```
0399
            FORMAT('<<CREAT MESSAGE>>IERR= ',14,' AT ',12)
            FORMAT('<<OPEN MESSAGE>>)[ERR= ',14,' AT ',12)
0600
                            MESSAGE>>IERR= '.I4.' AT '.I2)
0601
            FORMAT('<<LOCF
            FORMAT('<<CLOSE MESSAGE>>) IERR= ', 14, ' AT ', 12)
0602
            FORMAT('<<URITE MESSAGE>>) IERR= ', 14,' AT ', 12)
0603
            FORMAT('<<READF MESSAGE>>IERR= ',14,' AT ',12)
0604
0605
           FORMAT(' *****TIME OUT ERROR*****
           +/'******BUFFER OVERFLOW****', 16)
0606
0607
            END
0608
0609
            SUBROUTINE TAPE(TAPELU), TO CONTROL A TAPE LU BY PROGRAM
      C : AM A GEN PURPOSE ROUTINE TO DO HEAT THINGS TO ANY MAG-
0610
                                                                     C
0611
      C HETIC TAPE DRIVE.
0612
      C----
0613
            COMMON LU(5)
0614
            INTEGER TAPELU, IBUF(32)
9615
       50
            WRITE(LU,100) TAPELU
            FORNAT('----','ENTER AN OCTAL FUNCTION FOR LU ',12,' [04]'
0616
      100
0617
           +/' 0000 = CLEAR THE DEVICE'
           +/' 0100 = WRITE AN EDF'
0618
           +/' 0200 - BACKSPACE 1 RECORD'
0619
           +/' 0300 = FORWARD SPACE 1 RECORD'
0620
           +/' 0400 = REWIND'
0621
           +/' 0500 = WRITE AN INTER-RECORD GAP'
0622
           +/' 1300 = FORWARD 1 FILE'
0623
0624
           +/' 1400 = BACK 1 FILE'
           +/' 2700 = LOCATE FILE(CTU ONLY)'
0625
           +/' 7777 = READ & LIST A RECORD OF UP TO 32 WORDS'
0626
4627
           +/' 1111 = TO GET OUT OF THIS ROUTINE'/'-----')
0628
            READ(LU.101) IFC
0629
       101 FORMAT(04)
0630
            IF(IFC.NE.11118) GO TO 112
9631
                  WRITE(LU,113)
0632
       113
                  FORMAT( ' RETURNING TO MAIN PROGRAM')
0633
                  RETURN
           IF(IFC.EQ.2700B) GO TO 105
0634
       112
0635
            IF(IFC.EQ.77778) GO TO 110
0636
            I CODE=3
0637
            I 0P1 =0
           CALL EXEC(ICODE, TAPELU+IFC, IOP1)
0638
       103
0639
            GD TD 50
0640
            WRITE(LU, 106)
0641
            FORMAT(10X, 'ENTER FILE # TO LOCATE [12]')
0642
            READ(LU, +)[0P1
0643
            GO TO 103
0644
            CALL EXEC(1, TAPELU, IBUF, 32)
0645
            WRITE(LU, 111)(IBUF(I), I=1,32)
0646
            FORMAT(32A2)
0647
            GO TO 50
0648
            EHD
```

Ĺ

L

APPENDIX C - EXAMPLE: SIX on METHANOL PAN FIRE

REST1 T=00003 IS OH CROO041 USING 00056 BLKS R=0000

.

```
6 CM PAN FIRST DATA POINT
0001
                                         HORIZONTAL (M)
0002
                 VERTICAL (M)
                                                             SAMPLE INT. TIME (SEC)
0003
                 14526880
                                       -. 21625594E-02
                                                                .81100792E-01
                 31290054E-02
                                         43555414E-02
                                                                .46202880
3004
0005
                 73017776E-02
                                         57001272E-02
                                                                .38174719
9006
               - . 22369782E-02
                                        .40457380E-03
                                                                . 23592958
               -. 22369782E-02
0007
                                       - 71435757E-02
                                                                . 31948799
                 45054266E-02
                                       -. 21625594E-02
0009
                                                                 69795835
               -. 91599929E-03
                                         40457380E-03
                                                                . 42270720
3004
0010
                .63727051E-01
                                        .30249106E-02
                                                                .41615355
0011
                . 16927639
                                        .40457380E-03
                                                                .82329601E-01
0012
                . 16567457
                                        .17079869E-02
                                                                .10239999
0013
                .17667629E-02
                                        .40457380E-03
                                                                . 32767999
                . 3129v054E-02
0014
                                        .43555414E-02
                                                                46202880
                 31290054E-02
                                        .40457380E-03
                                                                 .78233600E-01
0015
0016
                .45054266E-02
                                       -. 21625594E-02
                                                                .71434236
                .41849620E-03
                                        .40457380E-03
                                                                . 57671678
0017
0018
                .31290054E-02
                                         17079869E-02
                                                                .14581758
0019
                .31290054E-02
                                        .17079869E-02
                                                                .11796479
0020
                 31290054E-02
                                         17079869E-02
                                                                .60293114
0021
                .31290054E-02
                                        .40457380E-03
                                                                .31293440
0022
                . 26740372
                                         40457380E-03
                                                                .38830078
0023
                .41849620E-03
                                       -. 21625594E-02
                                                                . 64225280
0024
                . 13424206
                                       -. 21625594E-02
                                                                . 24739838
0025
                .57635680E-01
                                       -. 88552607E-03
                                                                .31743996E-01
                                        .30249106E-02
                                                                . 44236797
0026
                . 43755323
0027
                .58962824E-02
                                        .43555414E-02
                                                                . 49479675
                                                                .53084159
0028
                 73017776E-02
                                        .30249106E-02
0029
               -. 22369782E-02
                                        .40457380E-03
                                                                . 17121279
                 44725734
                                         30249106E-02
0030
                                                                 50135040
0031
                 87221134E-02
                                        .57001272E-02
                                                                .41123837
0032
                10.004017
                                        .17079869E-02
                                                                .63487992E-01
0033
                9 1733093
                                        .30249106E-02
                                                                .10158080
0034
                4.6847963
                                        .40457380E-03
                                                                . 18022400
                14.246948
                                        .40457380E-03
                                                                .86425588E-01
0035
0036
                                       -. 88552607E-03
                13.454493
                                                                .10813439
0037
                .31290054E-02
                                        .43555414E-02
                                                                . 62586880
                                       -. 88552607E-03
                                                                .43212794E-01
0038
                6.5395241
0039
                4.7763882
                                        .43555414E-02
                                                                . 32767999
                6.6842194
0040
                                        .17079869E-02
                                                                . 32767999
0041
                7.6619654
                                         17079869E-02
                                                                .54886393E-01
                                        .30249106E-02
0042
                10.337299
                                                                .66764802E-01
0043
                .15864873
                                       -. 88552607E-03
                                                                49807358
0044
                .41849620E-03
                                        .30249106E-02
                                                                . 19496959
0045
                 31290054E-02
                                        .40457380E-03
                                                                . 83558393
0046
                .15522179
                                        .17079869E-02
                                                                .50135040
0047
                .55655666E-01
                                       -. 59169848E-02
                                                                . 28671998
0048
                .73017776E-02
                                        .30249106E-02
                                                                 50135040
0049
                .59640586E-01
                                        .30249106E-02
                                                                .70451200
0050
                .41849620E-03
                                       -. 21625594E-02
                                                                . 42270720
                 61670892E-01
0051
                                        .17079869E-02
                                                                 49807358
0052
                 41849620E-03
                                        .30249106E-02
                                                                .69468153
               -.91599929E-03
0053
                                       -. 88552607E-03
                                                                .65945596E-01
0054
                . 58962824E-02
                                        .30249106E-02
                                                                13926399
9055
                 99823773
                                        .40457380E-03
                                                                 31293440
0056
                . 15692824
                                       - .88552607E-03
                                                                . 55705595
0057
                .41849620E-03
                                        .40457380E-03
                                                                .66846716
0058
                . 45221126
                                        .40457380E-03
                                                                .80936956
```

```
0059
                 . 59608161
                                        -. 88552607E-03
                                                                  .12779519
0060
                 .17667629E-02
                                          .43555414E-02
                                                                  . 54722559
0061
                 . 12525600
                                          40457380E-03
                                                                  . 33587199
0062
                 31290054F-02
                                          30249106E-02
                                                                  .40468478
0063
                 .31290054E-02
                                          .57001272E-02
                                                                   38830078
0064
                -. 22369782E-02
                                          .30249106E-02
                                                                  . 27033597E-01
0045
                                         .40457380E-03
                  31290054E-02
                                                                   27361280
0066
                -. 73894039E-02
                                          17079869E-02
                                                                   11960319
0067
                -. 91599929E-03
                                          17079869E-02
                                                                   74055672
0068
                 . 44725734
                                         .17079869E-02
                                                                  .71434236
0069
                 .41849620E-03
                                         .30249106E-02
                                                                  .33095676
0070
                  58962824E-02
                                          43555414E-02
                                                                   24412158
0071
                  17667629E-02
                                          57001272E-02
                                                                   61276162
9972
                  17667629E-02
                                         .17079869E-02
                                                                  . 54722559
0073
                 . 37129261E-01
                                         .17079869E-02
                                                                  . 11468799
0074
                  45054266E-02
                                          57001272E-02
                                                                   64552951
0075
                  17667629E-02
                                         .17079869E-02
                                                                  . 29982716
0076
                 91599929E-03
                                          40457380E-03
                                                                   63242233
0077
                  45054266E-02
                                         .30249106E-02
                                                                  56688643
0078
                 .63727051E-01
                                         .43555414E-02
                                                                  .17940480
0079
                  41849620E-03
                                         .17079869E-02
                                                                 .72089601
0080
                 31290054E-02
                                          57001272E-02
                                                                  . 32604158
0081
                 17.288971
                                          88552607E-03
                                                                  . 14335999
0082
                 5.4457102
                                          30249106E-02
                                                                  . 14663678
0083
                 1.5851901
                                          30249106E-02
                                                                 .78643203E-01
0084
                 .45054266E-02
                                        -. 88552607E-03
                                                                  57343996
0085
                 .72924316
                                         .40457380E-03
                                                                 . 50790393
0086
                -. 91599929E-03
                                        - 21625594E-02
                                                                 . 31293440
0087
                 .16746795
                                          43555414E-02
                                                                 . 36864001
0088
                 . 61670892E-01
                                         .30249106E-02
                                                                 81592309
0089
                 .45054266E-02
                                         .17079869E-02
                                                                  57343996
0090
                - . 22369782E-02
                                         .30249106E-02
                                                                 .72089601
0091
                 .41849620E-03
                                         .40457380E-03
                                                                  86835191E-01
0092
                 .17667629E-02
                                         .17079869E-02
                                                                 .11468799
0093
                 .17667629E-02
                                         .40457380E-03
                                                                 . 21135360
0094
                 . 58962824E-02
                                         57001272E-02
                                                                 .71106553
0095
                 .16927639
                                         .70588645E-02
                                                                 . 37847036
0096
                 . 19425631
                                         .17079869E-02
                                                                  82903039
0097
                 . 41849620E-03
                                          30249106E-02
                                                                 .32358401E-01
0098
                 .17667629E-02
                                         .30249106E-02
                                                                  31129599
0099
                -.91599929E-03
                                        -.88552607E-03
                                                                  39157760
0100
                . 12235069
                                         .40457380E-03
                                                                 .66191351
0101
                 . 16567457
                                         .17079869E-02
                                                                 . 52428794
0102
                 .31290054E-02
                                         .30249106E-02
                                                                  66355199E-01
0103
                .31290054E-02
                                         .17079869E-02
                                                                 . 73072636
0104
                 .17667629E-02
                                         .17079869E-02
                                                                 . 28999680
0105
                  41849620E-03
                                         17079869E-02
                                                                  41779196
0106
                 41849620E-03
                                         .40457380E-03
                                                                 . 23756799
0107
                  45054266E-02
                                         .17079869E-02
                                                                 . 38993919
0108
                 41849620E-03
                                          30249106E-02
                                                                 .41287678
0109
                 .17667629E-02
                                         .17079869E-02
                                                                  83558393
0110
                .37641710
                                        -. 88552607E-03
                                                                 .42926079
0111
                .17110023
                                        -. 88552607E-03
                                                                  68812799
0112
               -. 22369782E-02
                                          17079869E-02
                                                                 .35717118
0113
                 .63727051E-01
                                         .30249106E-02
                                                                 .72089598E-01
0114
                . 35709697
                                         .40457380E-03
                                                                 . 12533760
0115
                . 13120064
                                          40457380E-03
                                                                 . 22282240
0116
                                         .40457380E-03
               -. 22369782E-02
                                                                 . 24248320
0117
                16567457
                                         .17079869E-02
                                                                 .16629758
0118
                .41949620E-03
                                         .17079869E-02
                                                                 69468153
```

•

1

1

•

.

4119	3:460033		
0119	. 36468273	. 30249106E-02	. 50462723
0120	82406759	88 5526 07E-03	12779519
9121	. 31290054E-02	.17079869E-02	. 15892479
0122			
5.5		21625594E-02	. 22609919
	22369782E-02	.40457380E-03	.43008000E-01
0124	. 59640586E-01	17079869E-02	. 82247674
0125	17479461	.43555414E-02	. 13107198
0126	19629845	.57001272E-02	. 70778871
0127	.31290054E-02	17079869E-02	. 40468478
0128	. 64159811	43555414E-02	25559038
0129	.31290054E-02		
	· —	.70588645E-02	. 73400319
0130	589 62824E -02	.70588645E-02	. 20479998
0131	40125716	.40457380E-03	. 56360960
0132	17667629E-02	.30249106E-02	.98303989E-01
0133			
		88552607E-03	. 19824639
0134	.45054266E-02	.70588645E-02	.46202880
0135	31290054E-02	.40457380E-03	. 44564480
0136	59640 586E-01	17079869E-02	. 27361280
		- 10751111E-01	. 73728001
133	41349629E-03	30249106E-02	.73728001
0139	.31290054E-02	. 57001272E-02	. 17694718
0140	59640586E-01	17079869E-02	64225280
	- 22369782E-02		
		30249106E-02	. 83558393
0142	62205803	.30249106E-02	. 15400958
0143	. 40558666	88552607E-03	. 78643191
0144	31290054E-02	.30249106E-02	. 44564480
0145	45054266E-02	17079869E-02	.70451200
0146	41849620E-03	.70588645E-02	. 43581438
0147	41849620E-03	.43555414E-02	.80936956
0148	.41849620E-03	.40457380E-03	. 47513598
0149	.31290054E-02		
		.30249106E-02	.93388796E-01
0150	.87221134E-02	88552607E-03	. 59310079
0151	31290054E-02	.43555414E-02	. 80936956
0152	.58962824E-02	.30249106E-02	.60293114
	91599929E-03	.40457380E-03	. 36208636
0154	17667629E-02	.40457380E-03	15646720
0155	.31290054E-02	. 17079869E-02	. 11304960
0156	41849620E-03 -	·.88552607E-03	. 39976960
0157	17667629E-02	.57001272E-02	36864001
0158	. 44725734	. 57001272E-02	. 54067194
0159	.73017776E-02	. 57001272E-02	. 71434236
0160	.31290054E-02	.43555414E-02	. 27033597E-01
0161	12235069	17079869E-02	. 13352960
0162	.31290054E-02	.30249106E-02	.77004790E-01
0163	65009533E-01	30249106E-02	. 50135040
0164	. 15522179	.17079869E-02	. 52101111
0165 -	.91599929E-03 -	. 88552607E-03	. 78315520
0166			
		·. 88552607E-03	. 28508157
0167	73017776E-02	.43555414E-02	. 73400319
0169	45054266E-02	.40457380E-03	. 77332473
0169 -	.91599929E-03 -	. 88552607E-03	. 10158080
	91599929E-03		
		.40457380E-03	. 82903039
0171	. 84765923	.40457380E-03	72499201E-01
0172	38448751	.70588645E-02	. 47513598
0173	41849620E-03	.40457380E-03	.95027193E-01
0174	.59640586E-01	.40457380E-03	. 45547515
0175	73017776E-02	. 43555414E-02	. 49151999
0176	.61670892E-01	.17079869E-02	.33382401E-01
0177	41849620E-03	.40457380E-03	. 26869756
0178	66202211	30249106E-02	. 56360960
· • •	· ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~ ~	OVETAINAT AT	. 54544784

```
0179
                 41849620E-03
                                       - 88552607E-03
                                                                31293440
0180
                 61670892E-01
                                        .40457380E-03
                                                                45875198
0181
                .91138351
                                         40457380E-03
                                                               .47513598
0192
                 31290054E-02
                                        .40457380E-03
                                                               . 33259517
0183
                 17110023
                                        .43555414E-02
                                                               . 24412158
0184
                .41849620E-03
                                        .43555414E-02
                                                               . 21954557
0185
                .45054266E-02
                                        .57001272E-02
                                                               . 66191351
0186
                 17667629E-02
                                         40457380E-03
                                                               . 16957438
0187
                91599929E-03
                                         40457380E-03
                                                               .53247996E-01
2188
                 41849620E-03
                                         17079869E-02
                                                                21299198E-01
0139
                 17667629E-02
                                       - 88552607E-03
                                                               .67502081
0190
                 45054266E-02
                                         40457380E-03
                                                               .80936956
0191
                .41849620E-03
                                        .43555414E-02
                                                               . 22609919
0192
               -. 61206026E-02
                                        .57001272E-02
                                                               . 34078717
0193
               -. 22369782E-02
                                        .30249106E-02
                                                               .14499840
0194
               -. 13548516E-01
                                        .30249106E-02
                                                               .69468153
0195
               -.91599929E-03
                                        .40457380E-03
                                                               . 39321595
0196
                . 28436261
                                        .30249106E-02
                                                               .64307198E-01
0197
                .45054266E-02
                                        .43555414E-02
                                                               . 33751041
9198
               -.91599929E-03
                                        .40457380E-03
                                                               .11960319
0199
                .17667629E-02
                                        .17079869E-02
                                                               . 23592958
0200
                . 59640586E-01
                                        .17079869E-02
                                                               .69140470
0201
                .41849620E-03
                                        .40457380E-03
                                                               .90112001E-01
0202
                . 27015656
                                       .30249106E-02
                                                               .78315520
0203
                .17667629E-02
                                      - 88552607E-03
                                                               . 12779519
0204
                . 98044619E-01
                                      -. 21625594E-02
                                                               .79953921
0205
                . 59640586E-01
                                       .17079869E-02
                                                               .79953921
0206
                .31290054E-02
                                       .30249106E-02
                                                               .72908789E-01
0207
0208
0209
0210
                 UNWEIGHTED MEAN VELOCITIES:
0211
                (CHTR 1)
                            .53307307 M/S
0212
                (CHTR 2>
                            .17951876E-02 H/S
0213
0214
                 RESULTANT VECTOR:
                 HAGNITUDE = 0.53307605
0213
0216
                 INCLUDED ANGLE = 0.19294965 FROM VERTICAL
0217
```

١,

D.

D

D

D

!..

Ĺ

19-84